

Products for Histology




Boekel Lighted Tissue Flotation Baths

Models 145951 and 145702

Features:

- Eliminate wrinkles and distortion during preparation of paraffin embedded tissue sections
- L.E.D. illuminated black background with transparent Borosilicate dish; allows specimen to be easily seen
- Digital microprocessor control with advanced safety over temperature protection
- Deluxe model includes Slide Dryer Block to assist in drying of slides and Tissue Orienter Block to smooth wrinkled sections



• Removable Borosilicate glass dish, Part Number: 904-0013

Accessories

• Borosilicate glass dish, Part Number: 904-0013



Boekel Slide Dispenser Model 1415:

- Stainless steel slide dispenser holds up to 54 3 in. x 1 in. slides
- Adjustable thickness
- Rubber cushion on bottom
- More specs on www.boekelsci.com



145951 Deluxe Lighted Tissue Flotation Bath



145702 Standard Lighted Tissue Flotation Bath

	Standard Lighted Tissue Flotation Bath	Deluxe Lighted Tissue Flotation Bath
Model Number	145702 (115 VAC) 145702-2 (230 VAC)	145951 (115 VAC) 145951-2 (230 VAC)
Overall Dimensions (Width x Length x Height)	36.2 cm x 36.2 cm x 11.4 cm (14.25 in x 14.25 in x 4.5 in)	36.2 cm x 36.2 cm x 11.4 cm (14.25 in x 14.25 in x 4.5 in)
Dryer Block	NA	5.1 cm x 12.7 cm (2 in x 5 in)
Orienter Block	NA	5.1 cm x 5.1 cm (2 in x 2 in)
Net Weight	4.7 kg (10.5 lb)	
Shipping Weight	5.4 kg (12 lb)	
Temperature Set Range	Ambient +5° C to 65° C	
Accuracy	+/- 1° C	
Uniformity at 55° C	+/- 1° C	
Dryer Block Temperature	NA	40° C +/- 3° C preset
Orienter Block Temperature	NA	55° C +/- 3° C fixed
Approvals	Certified to UL61010-1 and CSA C22.2	
Electrical	145702 (115 VAC 50/60 Hz 272 W) 145702-2 (230 VAC 50/60 Hz 272 W) (CE)	145951 (115 VAC 50/60 Hz 272 W) 145951-2 (230 VAC 50/60 Hz 272 W) (CE)

Boekel Standard Tissue Bath


Model 14792

Features:

- · Edge serves as drying surface for slides
- Surface is black porcelain enamel coated steel
- Equipped with On/Off switch, pilot light, thermostat and non-skid rubber feet
- · Rapid heat-up

Supplied With:

- · Bath Cover
- Thermometer Holder and Thermometer



Boekel Tissue Baths have a 2 year warranty on parts and workmanship

Boekel Slide Moat

Model 240000

Features:

- · 30 standard microscope slide capacity
- PID temperature controller for easy temperature selection
- Tempered glass lid which seals with integral gasket on heating base
- · Excellent temperature uniformity and stability



	Slide Moat Model	
Model Number	240000 (115 VAC), 240000-2 (230 VAC)	
Overall Dimensions (Width x Length x Height)	40.6cm x 34.9cm x 14.6cm (16" x 13 3/4" x 5 3/4")	
Net Weight	8.2 kg (18 lb)	
Shipping Weight	10.4 kg (23 lb)	
Capacity	30 - 1 x 3 inch slides (3 rows of 10)	
Temperature Range	Ambient + 5°C to 100°C	
Stability*	Typically +/- 0.2° C or better throughout operating temperatures	
Uniformity**	+/- 0.35° C at 37° C; +/- 0.7°C at 65° C	
Clearance between lid and slide pedestal	.48 in (1.0 cm)	
Heat-up Rate	Ambient to 65° C in 20 minutes	
Temperature Setting/Display	PID/LED	
Electrical	Model 240000: 115 VAC; 50/60 Hz; 200 W, (UL/ CSA)*** Model 240000-2: 230 VAC, 50/60 Hz, 155 W (CE)	



Boekel Slide Storage Cabinets

Model 141000

Features:

- Six cabinet drawers hold up to 5000 standard
 1 inch x 3 inch microscope slides
- · Cabinets can be interlocked up to 10 units high
- Plastic drawer dividers can be removed to store objects up to 5 cm (2 in) wide
- · Drawers slide smoothly and easily
- · Back of drawer acts as handle for transportation
- · Durable cold rolled steel with baked on epoxy finish



Accessories:

- Sponge block spacers
- Cabinet Base with Rubber Feet, Part Number: 141001
- Identification labels
 Drawer Only, Part Number: 141002

	Slide Storage Cabinet Model 141000	Cabinet Base Model 141001
Overall Dimensions (Width x Length x Height)	40 cm x 48.3 cm x 12.7 cm (15 3/4 in x 19 in x 5 in)	40 cm x 47.6 cm x 7 cm (15 3/4 in x 18 3/4 in x 2 3/4 in)
Net Weight	7.7 kg (17 lb)	1.3 kg (2 3/4 lb)
Shipping Weight	9.1 kg (20 lb)	1.8 kg (4 lb)

Note: Cabinets are stacked by utilizing built-in tabs in the bottom of each unit. Because of this, cabinets will "rock" and be unstable if placed on a bench-top without the use of a Cabinet Base, Part Number: 141001, underneath the first cabinet. See photo.

Boekel Wax Dispenser

Model 145600

Features:

- 5 Liter capacity, Approximately 1.25 Gallon
- · Heated valve ensures a consistent flow of molten wax
- Stainless steel outer casing
- Aluminum inner reservoir provides excellent thermal transfer and temperature control
- Dispenser valve allows precise control of flowing wax and can be locked in the open position if desired
- Inner reservoir screen helps to prevent large particulate from clogging the dispenser valve
- Rapid warm-up with an advanced safety overtemperature cut-off

	Wax Dispenser
Model Number	145600 (115 VAC), 145600-2 (230 VAC) (CE)
Overall Dimensions	21.6 cm x 40.6 cm (8 1/2 in x 16 in)
Temperature Range	40° C to 75° C +/- 2° C
Approvals	Certified to UL61010-1 and CSA C22.2
Electrical	115 VAC 50/60 Hz 250 W (CE)
Shipping Weight	6.3 kg (14 lb)



855 Pennsylvania Boulevard • Feasterville, PA 19053 USA 800-336-6929 • Phone: 215-396-8200 • Fax: 215-396-8264

E-mail: boekel-info@boekelsci.com • Web site: www.boekelsci.com

